



Your Vision Is Our Vision

As the Swiss leader in vision care, Ziemer develops premium eye care technology to deliver the best visual outcomes.



10 Million Eyes*

Ziemer technologies have supported more than 10 million eyes worldwide.



Innovative Pioneer


For more than 25 years, Ziemer has shaped the ophthalmic industry with innovations.



Global Presence

Present in more than 100 countries, Ziemer is trusted by leading professionals for precision solutions.

*Disclaimer: The figures referenced are derived from internal company data.



Unframe your life with AQUA PRK®



Scan QR Code

Learn more about the spa-like AQUA PRK treatment.

Ziemer Ophthalmic Systems AG
Allmendstrasse 11,
2562 Port, Switzerland

+41 32 332 70 50
marketing@ziemergroup.com

www.ziemereyecare.com/aquaprk

Wet Ablation

AQUA PRK

Glasses Off, Life On.

Experience the **AQUA** difference.

Issue Date 2026-04
Doc. No. AL4900-1402-01



Whether you're nearsighted or farsighted, with or without astigmatism, AQUA PRK® might be just the right treatment for you.

Spa-like Refractive Treatment



Comfort For Your Eyes

AQUA PRK® keeps your eyes naturally hydrated during surgery, designed to reduce dryness and enhance comfort compared to traditional PRK procedures.



Oderless and Quiet Treatment

Advanced technology creates an odorless and quiet environment, helping you feel at ease throughout your treatment.



Precise & Predictable Vision

Water-friendly wavelength aims for consistent results regardless of your corneal hydration level, so you can enjoy clear vision with confidence.

Your AQUA PRK® procedure

The goal of AQUA PRK surgery is to improve your vision by reshaping the cornea. Changing the cornea's shape adjusts its focusing power, helping your eye see clearly without the need for glasses or contact lenses.

01 | Personal Consultation & Eye Mapping

During your doctor's consultation, a personalized treatment plan is made after measuring and examining your eyes.

02 | Corneal Surface Removal

The epithelium is removed, typically using an alcohol solution, a mechanical instrument, or a specialized brush.

03 | Corneal Reshaping

A solid-state laser reshapes the cornea to improve the way light is focused, giving you clearer vision.

04 | Steady Recovery

A bandage-like soft contact lens is placed on the cornea to help protect the eye as it heals.

Could AQUA PRK® be right for you?

PRK is best suited for people who:

Enjoy contact/extreme sports, work in military, or have physical jobs

Prefer a procedure without cutting a corneal cap

Have a corneal structure that's better matched to PRK than LASIK

Advantages:

Safer for contact sports and high-impact activities

Free from corneal cap-related risks

Stable and predictable outcomes*

*Disclaimer: This statement pertains to all conventional PRK procedures

Consult your doctor or clinic staff to learn about any potential contraindications.

Disclaimer: The approved regulatory name of this module is PRK. Availability of this module may vary by country.